


# Tornier Perform™ Total Shoulder System

Properly size  
from your smallest to largest patients



## ◀ Stem length

Available in short,  
long and XL sizes

## Augmentation ▶

Defect-mimicking implants  
designed for bone preservation.



## ◀ Elevated by Blueprint

Surgeon-controlled  
3D planning software

stryker

# Improving your shoulder Arthroscopy

## Introducing **Tone Mode**



**Standard Mode +  
Tone mode Off**



**Standard Mode +  
Tone Mode On**

Designed to balance lighting across field of view by improving posterior lighting without compromising foreground detail.

## High Dynamic Range of **Suture Color**

- Helps provide clearer visualization under arthroscopy



Distinct anchor suture colors designed for visibility, complemented by HDR technology to improve contrast and detail.



# Tornier Perform

## Humeral System - Fracture

Do more for  
**fractures**



### Stem options

Distal stem diameter spans from  
**7mm to 17mm in 1mm increments**

### Bone graft window

Allows for tuberosity  
attachment to autograft



### Three fixation methods

- Cemented
- Cementless without screws
- Cementless with three screws

