



4.75 mm Corkscrew[®] FT Implant


Minimally Vented

- Double-loaded
 - 1.3 mm SutureTape with #2 tails
 - BioComposite screws
-
- Triple-loaded
 - #2 TigerTail[®] suture
 - BioComposite screws



Arthrex[®] 


When You Treat Shoulders...




*Knotless SutureTak®
Anchor*



*SwiveLock® Biceps
Tenodesis*



*Knotless PushLock®
Anchor*



*Superior Capsular
Reconstruction (SCR)*

...Think Arthrex

Arthrex ACP[®] Therapy Double Syringe System

Simple and Safe PRP Preparation


Benefits

- Autologous, regenerative therapy
- Quick preparation including 5 min. centrifugation at 1500 rpm
- Closed and sterile system for preparation and application

Mode of Action

- Inhibits potential inflammatory processes
- Initiates regeneration processes and pain reduction
- Supports the healing processes

Autologous
Conditioned
Plasma



www.arthrex.com

© Arthrex Inc. 2020. All rights reserved.

AD13-000285-en-US

Arthrex[®]