


ArthroFLEX[®] Dermal Allograft

Competitive Matrix






Arthrex[®] 

4.75 mm Corkscrew[®] FT Implant

Minimally Vented

- Double-loaded
- 1.3 mm SutureTape with #2 tails
- BioComposite screws

- Triple-loaded
- #2 TigerTail[®] suture
- BioComposite screws



Arthrex[®]

evBR13-003400-en-KR_B

Synergy^{ID}™


Near-Infrared Fluorescence 4K Imaging

See more than ever before with new image modes for customizable white light color augmentation, 4 times the number of pixels as HD, a multi-sensory camera design, 10-bit precision color reproduction, and direct LED integration. Using endoscopic imaging, the Synergy^{ID}™ UHD 4K system offers dynamic range and depth of light for unrivaled contrast and depth of field for optimal sharpness.


At the touch of a button, new visualization options include standard visible light imaging with near-infrared (NIR) overlay, grayscale visible light imaging with NIR overlay, or NIR imaging-only modes. This valuable modularity provides immediate access to critical features with flexibility to scale investment in the system to meet future needs.

Reveal the Power of Fluorescence


Experience the Synergy^{ID}™ all-in-one console's* outstanding 4K video with advanced visualization, fluorescence imaging, LED lighting, image management, OR integration, and an easy-to-use tablet interface that allows staff to focus on patients instead of equipment. Multiple modes and color options put premier customized visualization at your fingertips.



Standard visible light



Standard visible light
+ NIR overlay



Grayscale visible light
+ NIR overlay



NIR only



4K
Synergy ID™

Arthrex®
[Synergy.Arthrex.com](https://www.Synergy.Arthrex.com)

© 2021 Arthrex, Inc. All rights reserved. AD13-000304-en-US_A